

Kayıtlı

- Yüksek Doğruluk
- Kullanıcı Dostu
- Lis'e Bağlanabilme
- Kolay Kullanım
- Hızlı Sonuç
- Dahili Yazıcı ile Rapor Alabilme

Fluorescence Immunoassay Analyzers

Tek Kanallı



OPENLAB-600

Portatif



OPENLAB-301



OPENLAB-302

6 Kanallı



OPENLAB-400

12 Kanallı



OPENLAB-800

İmmunoflorasan Testleri



İMMUNOFLORESAN CİHAZLARI



İmmunofloresan Cihazları İmmunofloresan Testler

TEST MENÜSÜ


Turuncu ile işaretlenmiş ürünlerde: Parmak ucu tam kan numune alımı mümkündür.

Kat. No.	Ürün	Numune	Cut-off	Ölçüm Aralığı	Ambalaj	CE Durumu	Klinik Performans	Okuma Süresi
Enflamasyon Belirteçİ								
OP-CRP-402	CRP	Tam Kan Serum Plazma	hsCRP: 1 mg/L, CRP: 10 mg/L	0.5-200 mg/L	10T/25T	CE	R= 0.986	3 Dak.
OP-PCT-402	PCT	Tam Kan Serum Plazma	0.5 ng/mL	0.1-50 ng/mL	10T/25T	CE	R= 0.981	15 Dak.
OP-IL6-402	IL-6	Tam Kan Serum Plazma	10 pg/mL	3.0-5000 pg/mL	10T/25T	CE	R= 0.996	15 Dak.
OP-SAA-402	SAA	Tam Kan Serum Plazma	10 mg/L	1-300 mg/L	10T/25T	CE	R= 0.990	15 Dak.
Tümör Belirteçİ								
OP-PSA-302	PSA	Kan Plazma	4 ng/mL	2-100 ng/mL	10T/25T	—	R= 0.990	15 Dak.
OP-PSA-402	PSA	Tam Kan Serum Plazma	4 ng/mL	2-100 ng/mL	10T/25T	—		
OP-FPSA-302	Free PSA	Kan Plazma	/	0.1-30 ng/mL	10T/25T	—	R= 0.994	15 Dak.
OP-AFP-302	AFP	Kan Plazma	20 ng/mL	5-400 ng/mL	10T/25T	CE	R= 0.981	15 Dak.
OP-CEA-302	CEA	Kan Plazma	4 ng/mL	1-500 ng/mL	10T/25T	CE	R= 0.976	15 Dak.
OP-FOB-602	FOB	Gaita	50 ng/mL	10-1000 ng/mL	10T/25T	CE	R= 0.992	15 Dak.
OP-CA211-402	CYFRA21-1	Tam Kan Serum Plazma	<3.3 ng/mL	1-300ng/mL	10T/25T	Non-CE	R=0.990	15 Dak.
Kardiyak Belirteçİ								
OP-CTI-402	cTnI	Tam Kan Serum Plazma	0.5 ng/mL	0.1-40 ng/mL	10T/25T	CE	R=0.989	15 Dak.
OP-CTNT-402	cTnT	Tam Kan Serum Plazma	0.5 ng/mL	0.2-40 ng/mL	10T/25T	CE	R=0.9923	15 Dak.
OP-BBNP-402	NT-proBNP	Tam Kan Serum Plazma	0.45 ng/mL	0.3-22 ng/mL	10T/25T	CE	R=0.990	15 Dak.
OP-CKMB-402	CK-MB	Tam Kan Serum Plazma	5 ng/mL	0.2-75 ng/mL	10T/25T	CE	R=0.989	15 Dak.
OP-HFA-402	H-FABP	Tam Kan Serum Plazma	8 ng/mL	1-120 ng/mL	10T/25T	CE	R=0.989	15 Dak.
OP-MYO-402	Myoglobin	Tam Kan Serum Plazma	90 ng/mL	5-200 ng/mL	10T/25T	CE	R=0.990	15 Dak.
OP-AMI-435	Troponin I/Myoglobin/CK-MB (3 in 1)	Tam Kan Serum Plazma	cTnI: 0.5 ng/mL MYO: 90 ng/mL CK-MB: 5 ng/mL	0.1-40 ng/mL 5-200 ng/mL 0.2-75 ng/mL	10T/25T	CE	cTnI: R=0.997 MYO: R=0.993 CK-MB: R=0.997	15 Dak.
OP-AMIT-435	Troponin T/Myoglobin/CK-MB Combo (3 in 1)	Tam Kan Serum Plazma	cTnT:0.2 ng/mL MYO:90 ng/mL CK-MB: 5 ng/mL	1-40 ng/mL 5-200 ng/mL 0.2-75 ng/mL	10T/25T	CE	cTnT: R=0.988 MYO:R=0.980 CK-MB: R=0.986	15 Dak.
OP-CPC-435	Troponin I/NT-proBNP/ D-Dimer Combo (3 in 1)/Serum Plazma	Tam Kan	cTnI: 0.5 ng/mL D-Dimer: 0.5mg/L NT-proBNP: 0.45 ng/mL	0.1-40 ng/mL 0.1-10 mg/L 0.3-22 ng/mL	10T/25T	CE	cTnI: R=0.985 D-Dimer: R=0.988 NT-proBNP: R=0.991	15 Dak.
Koagülasyon								
OP-DDM-402	D-Dimer	Tam Kan/Plazma	0.5 mg/L	0.5 mg/L	10T/25T	CE	R= 0.991	15 Dak.
Endokrinoloji								
OP-FSH-402	FSH	Tam Kan Serum Plazma	/	1-100 mIU/mL	10T/25T	CE	R= 0.987	15 Dak.
OP-LH-402	LH	Tam Kan Serum Plazma	/	1-300 mIU/mL	10T/25T	CE	R= 0.988	15 Dak.
OP-T3-302	T3	Kan Plazma	1.23-3.08 nmol/L	0.62-9.24 nmol/L	10T/25T	CE	R= 0.975	15 Dak.
OP-T3-402	T3	Tam Kan Serum Plazma	1.23-3.08nmol/L	0.62-9.24 nmol/L	10T/25T	CE	R=0.975	15 Dak.
OP-FT3-302	FT3	Kan Plazma	2.6-5.7 pmol/L	1.5-46 pmol/L	10T/25T	CE	R= 0.996	15 Dak.
OP-FT3-402	FT3	Tam Kan Serum Plazma	2.6-5.7 pmol/L	1.5-46 pmol/L	10T/25T	CE	RR=0.996	15 Dak.
OP-T4-302	T4	Kan Plazma	65-155 nmol/L	12.87-300 nmol/L	10R/25T	CE	R=0.984	10 Dak.
OP-T4-402	T4	Tam Kan Serum Plazma	65-155nmol/L	12.87-300 nmol/L	10T/25T	CE	R=0.976	15 Dak.
OP-FT4-302	FT4	Kan Plazma	6.7-16.4 pmol/L	1-100 pmol/L	10T/25T	CE	R= 0.994	15 Dak.
OP-FT4-402	FT4	Tam Kan Serum Plazma	6.7-16.4pmol/L	1-100 pmol/L	10T/25T	CE	R=0.989	15 Dak.
OP-TES-302	Testosterone	Kan Plazma	Female:0-2 ng/mL, Male:2.6-12 ng/mL	1.0-20 ng/mL	10T/25T	CE	R= 0.975	15 Dak.
OP-BHCG-302		Kan Plazma	25 mIU/mL	5-10,000 mIU/mL	10T/25T	CE	R= 0.994	15 Dak.
OP-BHCG-402		Tam Kan Serum Plazma	5 mIU/mL	2-200,000 mIU/mL	10T/25T	CE	R= 0.995	15 Dak.
OP-P4-302	Progesterone(P4)	Kan Plazma	/	1.4-60 ng/mL	10T/25T	CE	R=0.975	15 Dak.
OP-FER-302	Ferritin	Kan Plazma	Male:30-400 ng/mL, Female:13-150 ng/mL	1.0-500.0 ng/mL	10T/25T	CE	R=0.984	15 Dak.
OP-TSH-302	TSH	Kan Plazma	/	1.0-200 ng/mL	10T/25T	CE	R= 0.984	15 Dak.
OP-PRL-402	PRL	Tam Kan Serum Plazma	Female (Non-pregnant): 2-25 ng/mL Female (Pregnant): 10-209 ng/mL Male:2 - 18 ng/mL	1.0 - 200 ng/mL	10T/25T	CE	R= 0.951	15 Dak.
OP-AMH-402	AMH	Tam Kan Serum Plazma	/	0.05-25 ng/mL	10T/25T	CE	R= 0.972	15 Dak.
OP-COR-302	Cortisol	Kan Plazma	140 - 690 nmol/L	50 nmol/L - 1000 nmol/L	10T/25T	CE	R= 0.981	15 Dak.
Netroloji								
OP-ALB-102	mAlb	İdrar	20mg/L	5-300 mg/L	10T/25T	CE	R= 0.990	10 Dak.
FL-CYSC-402	CysC	Tam Kan Serum Plazma	1.1 mg/L	0.1-9.0 mg/L	10T/25T	CE	R= 0.990	10 Dak.
OP-BMG-402		Tam Kan Serum Plazma	2.2 mg/L	0.2-20 mg/L	10T/25T	CE	R=0.990	10 Dak.
OP-NGAL-102	N-GAL	İdrar	10 ng/mL	10-1500 ng/mL	10T/25T	CE	R=0.992	10 Dak.
Diabet								
OP-HBA-402	HbA1c	Tam Kan	4-6%	4-14.5%	10T/25T	CE	R= 0.981	10 Dak.
Otoimmün								
OP-IGE-402	İgE	Tam Kan Serum Plazma	100 IU/mL	20-1600 IU/mL	10T/25T	CE	R= 0.983	15 Dak.
OP-IGA-402	İgA	Tam Kan Serum Plazma	0.7-4.3 g/L	0.1-7.0g/L	10T/25T	CE	R=0.987	15 Dak.
OP-RFM-402	Rheumatoid Factor (RF) IgM	Tam Kan Serum Plazma	20IU/mL	8-200 IU/mL	10T/25T	CE	R= 0.969	15 Dak.

! Clinical Performans, Kararlılık, Koriyasyon Katsayısı, Sensitivite, Spesifite (Doğruluk)

Kat. No.	Ürün	Numune	Cut-off	Ölçüm Aralığı	Ambalaj	CE Durumu	Klinik Performans	Okuma Süresi
Enfeksiyon Hastalıkları								
OP-STA-502	Strep A	Swab	/	Qualitative	10T/20T	CE	Sensitivity:97.8%/Specificity:100%/Accuracy:98.7%	15 Dak.
OP-FLU-502	FLU A+B	Swab	/	Qualitative	10T/20T	CE	Nasopharyngeal Swab: FLU A: Sensitivity:98.2%/Specificity:99.0%/Accuracy:98.7% FLU B: Sensitivity:97.8%/Specificity:98.9%/Accuracy:98.6% Throat Swab: FLU A: Sensitivity:97.6%/Specificity:99.3%/Accuracy:98.5% FLU B: Sensitivity:94.9%/Specificity:99.5%/Accuracy:98.2% Nasal Aspirate: FLU A: Sensitivity:100.0%/Specificity:99.2%/Accuracy:99.3% FLU B: Sensitivity:97.8%/Specificity:99.4%/Accuracy:98.8%	15 Dak.
OP-RSV-502	RSV	Swab	/	Qualitative	10T/20T	CE	Sensitivity:97.8%/Specificity:100%/Accuracy:98.7%	15 Dak.
OP-NCP-502	COVID-19 Antijen	Swab	/	Qualitative	20T	CE	Sensitivity:96.43%/Specificity:99.29%/Accuracy:97.53%	10 Dak.
OP-DEN-402	Dengue IgG/IgM	Tam Kan Serum Plazma	/	Qualitative	10T/25T	CE	Dengue IgG Sensitivity:91.07%/Specificity:92.99%/Accuracy:92.49% Dengue IgG Sensitivity:89.47%/Specificity:95.51%/Accuracy:92.49%	15 Dak.
OP-DES-402	Dengue NS1	Tam Kan Serum Plazma	/	Qualitative	10T/25T	CE	Sensitivity:96.8%/Specificity:91.9%/Accuracy:93.4%	15 Dak.
OP-NCP-402	COVID-19 IgG/IgM	Tam Kan Serum Plazma	/	Qualitative	10T/25T	CE	Sensitivity:96.7%/Specificity:96.9%/Accuracy:96.9%	15 Dak.
OP-CSG-402	COVID-19 S-RBD IgG Antikor	Tam Kan Serum Plazma	/	Semi-quantitative	10T/25T	CE	Sensitivity:97.72%/Specificity:99.25%/Accuracy:98.31%	15 Dak.
OP-CIC-525	COVID-19 Ag/Flu A+B Combo	Nazal Swab	/	Qualitative	20T	CE	COVID-19: Sensitivity:97.14%/Specificity:99.05%/Accuracy:98.58% FLU A: Sensitivity:97.87%/Specificity:98.71%/Accuracy:98.39% FLU B: Sensitivity:98.80%/Specificity:98.83%/Accuracy:98.82%	15 Dak.
OP-SYP-402	Syphilis Antikor	Tam Kan Serum Plazma	/	Qualitative	10T/25T	CE	Sensitivity:>99.9%/Specificity:99.3%/Accuracy:99.6%	15 Dak.
OP-SP-102	Streptococcus pneumoniae	İdrar	/	Qualitative	10T/25T	CE	Sensitivity:93.3%/Specificity:98.9%/Accuracy:98.1%	15 Dak.
OP-LEG-102	Legionella pneumophila Antijen	İdrar	/	Qualitative	10T/25T	CE	Sensitivity:97.1%/Specificity:>99.9%/Accuracy:99.0%	15 Dak.
OP-CSG-402B	COVID-19 S-RBD IgG Quantitative	Tam Kan Serum Plazma	/	Qualitative	10T/25T	CE	Sensitivity:97.6%/Specificity:99.0%/Accuracy:98.8%	15 Dak.
OP-HP-602	H. pylori Antijen	Gaita	< 2.5 ng/mL	2-1500 ng/mL	10T/25T	CE	Sensitivity:98.00%/Specificity:99.00%/Accuracy:97.00%	15 Dak.
OP-MAL-435	Malaria Pf./P.v./Pan /Antijen	Tam Kan	/	Qualitative	10T/25T	CE	Sensitivity:98.7%/Specificity:99.3%/Accuracy:99.1%	15 Dak.
OP-TYP-402	Typhoid Antikor	Tam Kan Serum Plazma	/	Qualitative	10T/25T	CE	IgM: Sensitivity:94.3%/Specificity:99.3%/Accuracy:98.8% IgG: Sensitivity:90.0%/Specificity:99.7%/Accuracy:99.0%	15 Dak.
OP-TET-402	Tetanus Antikor	Tam Kan Serum Plazma	0.1 IU/mL	0.05-10 IU/mL	10T/25T	CE	/	15 Dak.
OP-ROT-602	Rotavirus Antijen	Gaita	/	Qualitative	10T/25T	—	Sensitivity:98.00%/Specificity:96.00%/Accuracy:97.00%	15 Dak.
OP-HBSG-402	HBsAg Antijen	Tam Kan Serum Plazma	<0.5IU/mL	0.5-200 IU/mL	10T/25T	—	R= 0.991	/
OP-HBSB-402	HBsAb Antikor	Tam Kan Serum Plazma	<10mIU/mL	5-1000 mIU/mL	10T/25T	—	R= 0.981	15 Dak.
OP-HCV-402	HCV Antikor	Tam Kan Serum Plazma	/	Qualitative	10T/25T	—	Sensitivity:>99.9%/Specificity:99.00%/Accuracy:99.42%	15 Dak.
OP-MXG-402	Monkeypox Virus Antijen	/	/	Qualitative	20T	CE	/	15 Dak.
OP-HIV-402	HIV1/2 Antikor	Tam Kan Serum Plazma	/	Qualitative	10T/25T	—	Sensitivity:96.8%/Specificity:97.8%/Accuracy:97.5%	15 Dak.
OP-HIG-402	HIV P24 Antijen	Tam Kan Serum Plazma	/	Qualitative	10T/25T	—	/	15 Dak.
OP-ASO-302	Anti-streptolysin O (ASO)	Kan Plazma	Adults:<200 IU/mL Children:<250 IU/mL	25-800 IU/mL	10T/25T	CE	R= 0.972	15 Dak.
Diğerleri								
OP-CAL-602	Calprotectin	Gaita	/	Qualitative	10T/25T	CE	R= 0.971	15 Dak.
OP-VD-302	Vitamin D	Kan Plazma	>30 ng/mL	5-100 ng/mL	10T/25T	CE	R=0.978	15 Dak.
Uyuşturu Paneli								
OP-MET-102	MET	İdrar	1000 ng/mL	100-5000 ng/mL	10T/25T	CE	/	5 Dak.
OP-MOP-102	MOP	İdrar	/	Quantitative	10T/25T	CE	/	5 Dak.
OP-COC-102	COC	İdrar	300 ng/mL	50-3000 ng/mL	10T/25T	CE	/	5 Dak.
OP-AMP-102	AMP	İdrar	1000 ng/mL	300-3000 ng/mL	10T/25T	CE	/	5 Dak.
OP-BZO-102	BZO	İdrar	300 ng/mL	150-1500 ng/mL	10T/25T	CE	/	5 Dak.
OP-APVP-102	APVP	İdrar	/	Quantitative	10T/25T	CE	/	5 Dak.
OP-BAR-102	BAR	İdrar	300 ng/mL	50-1200 ng/mL	10T/25T	CE	/	5 Dak.
OP-THC-102	THC	İdrar	50 ng/mL	25-150 ng/mL	10T/25T	CE	/	5 Dak.
OP-MDMA-102	MDMA	İdrar	/	Quantitative	10T/25T	CE	/	5 Dak.
OP-MTD-102	MTD	İdrar	300 ng/mL	150-1500 ng/mL	10T/25T	CE	/	5 Dak.
OP-PCP-102	PCP	İdrar	/	Quantitative	10T/25T	CE	/	5 Dak.
OP-OPI-102	OPI	İdrar	2000 ng/mL	300-6000 ng/mL	10T/25T	CE	/	5 Dak.
OP-DOA-125-195	Multi-Drug 2-9 Drugs	İdrar	THC/COC/MET/ BAR/OPI/PCP/MDMA	See Insert	10T/25T	CE	/	5 Dak.

Klinik örneklerin başka bir marka ticari test kiti ile karşılaştırıldığında korelasyon katsayısı.



Uyuşturu Testi (THC) Test Raporu

Enfeksiyon Hastalıkları (Dengue) Test Raporu